Ammonium Decavanadate-Synthesis and Analysis 1 TSUD Department of Chemistry

In this experiment, we will synthesize ammonium decavanadate, $(NH_4)_6V_{10}O_{28}\cdot 6H_2O$, and then confirm the vanadium content by titration with permanganate.

Synthesis of Ammonium Decavanadate

Weigh 3.0 g of ammonium metavanadate (NH_4VO_3) into a 250 mL Erlenmeyer flask. Add 100 mL of distilled water and a stir bar and heat to just below boiling, while stirring constantly, until most of the solid has dissolved. Filter the solution and add, with stirring, 4 mL of 50% aqueous acetic acid. Add 150 ml of 95% aq ethanol and cool in an ice bath for 15 minutes. Suction filter the orange product and wash twice with 15 mL portions of ice-cold ethanol. Continue suction for 5 minutes and then transfer the product to a tared beaker. Cover loosely and store for at least two days to allow the crystals to air dry. Weigh after drying to determine the yield.

Vanadium Analysis

The V in the metavanadate will be reduced to V^{4+} and then titrated with MnO_4^{1-} according to the following equations:

$$V_{10}O_{28}^{6-} + 5H_2SO_3 + 26H^{1+} = 10VO^{2+} + 5H_2SO_4 + 13H_2O$$

 $5VO^{2+} + MnO_4^{1-} + H_2O = 5VO_2^{+} + Mn^{2+} + 2H^{1+}$

Weigh two 0.3xxx g samples of your product into 250 mL flasks and dissolve each in 40 mL of 1.5 M H_2SO_4 (warm if necessary). In the fume hood, add 50 mL of distilled water and 1 g of NaHSO₃. Swirl gently to dissolve and then allow to stand for 5 minutes. After standing, heat to a gentle boil and maintain for 15 minutes to remove excess H_2SO_3 . Titrate the solutions with 0.02xxx M MnO_4^{1-} while still warm. The equivalence point is the first persistent darkening of the yellow solution.

¹ adapted from D. C. Harris, Exploring Chemical Analysis, 2e, Freeman, 2001, pp 507-509.